

Laboratory Refrigerator

ตู้เย็นสำหรับห้องปฏิบัติการ 2-8°C

P701

Temperature Control

- Programable Controller
- Forced-air circulation system ensuring constant cabinet temperature
- Temperature setting range can selected by user from 2-8 °C

Ergonomic Design

- Multilevel adjustable shelves to keep pharmacy of various dimensions.
- LED display for easy observation.
- CFC-free insulator and refrigerator.
- Inner light.
- 4 castors easy to move.

Alarm and safety function

- Audible alarms.
- High and low temperature alarm.



Technical Data

Specification

Model P701

Capacity	700 Liter (24 cu.ft.)
Temperature Control	Programmable Controller
Temperature Range	2-8°C with 0.1°C resolution
Ambient Temperature	Not over 30 °C
Display	LED
Material internal and external	Powder paint
Alarm	Audible and Visual alarms Upper and Lower alarms Error at control sensor alarm
Type Door	2 glass doors
Inner light	LED x1 Set
Material Shelf	5 x 2 with PE Coating
Defrost	No frost
Refrigerant	R134a (CFC-Free)
Dimension (W x D x H)	Internal 1020 x 500 x 1380 mm
	External 1100 x 600 x 2000 mm.
Weight	115 Kg.
Power	220-240 Volt/50Hz
Warranty	1 Year